



Tedis Mobile App User Guide

Patient

Generated on: 2026-08-01


Contents

- | | |
|--------------------------------|----------------------|
| 1. Today | 8. Devices |
| 2. Observations | 9. Pending Alerts |
| 3. Observation Form | 10. Chats |
| 4. Image Source Selection | 11. Profile Details |
| 5. Warfarin | 12. Settings |
| 6. Weekly Dose Schedule Detail | 13. Report a Problem |
| 7. Medication Schedule | 14. About |

STEP 1

Today


The main screen you see when you open the app. The top card shows your latest INR value, the measurement date, and your next control date. Below it, today's medications are listed (dose, when and how to take it, with the doctor's note if any); tap "Taken" once you've taken a medication — this record lets your doctor track your medication adherence.



STEP 2

Observations

Shows the history of every measurement you've recorded (INR, weight, blood pressure) in chronological order. Use the "Add Observation" button in the bottom-right to enter a new measurement, or tap a record to see its detail. The data here is exactly what your doctor sees.



STEP 3

Observation Form

The form for entering a new measurement; fields shown adapt to your device configuration (INR, blood pressure, weight, height, or VAD device parameters). "Add Image" lets you attach a driveline (cable exit site) photo along with the observation — this helps your doctor remotely monitor the wound site.

 Observation Form

FLOW RATE

WATT

RPM

PI

INR

SYSTOLIC PRESSURE

DIASTOLIC PRESSURE

WEIGHT

HEIGHT


 Add Image

 Save

STEP 4

Image Source Selection

Opened by the "Add Image" button; lets you choose between taking a photo directly with the camera or picking one from your gallery. The app warns you beforehand about photos that may be too blurry to use.



< Observation Form

FLOW RATE

WATT

RPM

PI

INR

SYSTOLIC PRESSURE

DIASTOLIC PRESSURE

WEIGHT

HEIGHT

+ Add Image

Camera

Gallery

STEP 5

Warfarin

Shows the warfarin plan your doctor has defined for you: medication name, weekly dose distribution, your target INR range, and your next control date. Tap the weekly dose strip to see exactly how much to take each day (including any days marked "skip"). Remember: your dose is always determined by your doctor — the app only reminds you.

☰

Warfarin Treatment

<

August 2026

2 Weeks

>

Mon	Tue	Wed	Thu	Fri	Sat	Sun
27	28	29	30	31	1	2
3	4	5	6	7	8	9

Treatment List

Take in the evening, after a meal.
Doctor's note: Test ilaç planlama

VARFARİN 5 MG TABLET, 50
ADET
5.0 mg (1 full tablet)

Taken

Weekly Dose Schedule

Weekly Total: 35 mg

Mon	Tue	Wed	Thu	Fri	Sat	Sun
5	5	5	5	5	5	5

Today

Observations

Warfarin

Chats

STEP 6

Weekly Dose Schedule Detail

Opened by tapping the weekly dose strip on the Warfarin screen; shows the planned dose (in mg and tablet count) for each day of the week, plus the weekly total. Today is highlighted, and days with no dose are labeled "Skip".

Warfarin Treatment

< August 2026 2 Weeks >

Mon

Tue

Wed

Thu

Fri

Sat

Sun

27

28

29

30

31

1

2

3

4

5

6

7

8

9

Treatment List

i

Take in the evening, after a meal.
Doctor's note: Test ilaç planlama

VARFARİN 5 MG TABLET, 50
ADET
5.0 mg (1 full tablet)

Taken

Weekly Dose Schedule

Monday	5 mg (1 full tablet)
Tuesday	5 mg (1 full tablet)
Wednesday	5 mg (1 full tablet)
Thursday	5 mg (1 full tablet)
Friday	5 mg (1 full tablet)
Saturday	5 mg (1 full tablet)
Sunday	5 mg (1 full tablet)
Weekly Total: 35 mg	

STEP 7

Medication Schedule

Shows the plan for other medications you take regularly besides warfarin (blood pressure, heart failure medications, etc.) and your past "taken" records.



Medication Schedule



August 2026

2 Weeks



Mon	Tue	Wed	Thu	Fri	Sat	Sun
27	28	29	30	31	1	2
3	4	5	6	7	8	9

Treatment List

Today

Observations

Warfarin

Chats

STEP 8

Devices

Shows information about your implanted VAD device (type, serial number, location, registration date) — this is entered by your doctor; you can only view it.

☰

Devices

Right

Active

 HEARTMATE III
HEARTMATE 3

#

MLP_9999
Serial No

 12/7/2026
Device Installation Time

 20 Day(s)
Device's Total Usage Time

 Test amaçlı ikinci cihaz ve cihaz
kurumu seçimi
Description

Left

Inactive

 HEARTWARE
HEARTWARE

#

MLP-00000
Serial No

 Today

 Observations

 Warfarin

 Chats


Tedis — Heartsapiens Medikal Teknoloji A.Ş.

10

STEP 9

Pending Alerts


Lists clinical alerts about you that haven't been processed yet (e.g. an observation value outside the target range). If there are none, you'll see a green checkmark with "No pending alerts" — meaning there's currently nothing open for your doctor to follow up on.



STEP 10

Chats

The screen where you can message your doctor or institution directly. This channel is meant for short, non-urgent questions and updates — in an emergency, always contact your healthcare provider directly.



STEP 11

Profile Details

Shows your account information (first name, last name, institution, email, phone number) as read-only. If any of this needs to change, contact your institution.

 Profile Details

Name
App Test

Surname
User

Organisation
TEST KURUMU

E-mail
tedis@heartsapiens.art

Phone
+90(555)555555

 Today

 Observations


 Warfarin

 Chats

STEP 12

Settings

Toggle push notifications on/off, choose the app language (8 languages supported) and theme (Light/Dark/System). The app version number is shown at the bottom.



STEP 13

Report a Problem

A form for reporting an issue directly to the development team. You can write a subject and details, and optionally attach a screenshot.

☰

Report a Problem

Subject

Screenshot (Optional)

📷

Report Detail

Report

Today

Observations

Warfarin

Chats

STEP 14

About

Shows the app's version information, the developer company, and support/contact channels (email, website).

