



Tedis Web Application User Guide

Patient

Generated on: 2026-08-01

Contents

1. Dashboard
2. Pending Alerts
3. Patient Detail — General Information
4. Observation Values Tab
5. Notifications Tab
6. My Account

STEP 1

Dashboard


The Dashboard is the first screen you see after logging in to the Tedis web application. For patient accounts this screen does not contain any personal health data; it only shows a technical summary of overall system health (average CPU/HDD/RAM usage). From the left-hand menu you can reach "Pending Alerts" and your "Patient Detail" page. Clicking your name in the top-right corner opens your account settings, the flag icon lets you switch language, and the exit icon on the far right logs you out safely.



STEP 2

Pending Alerts

This screen lists clinical alerts about you that have not yet been processed — for example an observation value (such as INR) outside the target range, or a missed warfarin dose notification. If there are no alerts, you will see a green checkmark with "No pending alerts", meaning there is currently nothing open that your care team needs to act on regarding you. If an alert exists, clicking the row shows its type, the related observation, and its current status. This list is for information only — marking an alert as processed is done by your doctor or care team, not by you.



STEP 3

Patient Detail — General Information

The "Patient Detail" link in the left menu takes you to your own patient record. The top section shows your institution, ID number, name/surname, email, phone, date of birth, gender, nationality and target INR range in read-only form — only your doctor or an authorized staff member at your institution can change these; you cannot edit them. The "Allow Changes" lock is disabled for patient accounts. Options such as "Pause/Resume Notifications" in the top-right are specific to the doctor/manager roles and do not appear on the patient screen. Below, only two tabs relevant to you are shown: Observation Values and Notifications.

**Tedís**
Treatment Monitoring System

Version: 2026.08.01

[Dashboard](#)

[Pending Alerts](#)

[Patient Detail](#)

[App Test](#) EN [↗](#)

[Back](#)

Patient App Test User [Allow Changes](#)

Institution

TEST KURUMU

ID Number

12345678900

First Name

App Test

Last Name

User

Email

tedis@heartsapiens.art

Phone

+90(555)555555

Date of Birth

01.01.2000

Gender

☐ Female ☒ Male ☐ Other

Nationality

TR - Türkiye

Target INR Range

☒ 2.00-3.00 ☐ 2.50-3.50

Description

[Observation Values](#)

Tedís © 2026 Treatment Monitoring System

[Privacy Policy](#)

STEP 4

Observation Values Tab

This tab lists every observation value recorded for you through the mobile app or your devices (such as INR, weight, blood pressure, or VAD device parameters). You can search past records using the Start/End Time, Observation Type, and value range (Greater/Less/Equal/Between) fields. This screen is read-only; to add a new observation value you need to use the mobile app. The web interface is designed for reviewing your history and seeing the same data your doctor sees.

Tedis
Treatment Monitoring System
Version: 2026.08.01

[Dashboard](#)
[Pending Alerts](#)
[Patient Detail](#)

[App Test](#) [EN](#) [↗](#)

[Back](#)

Patient App Test User [Allow Changes](#)

Institution
TEST KURUMU

ID Number
12345678900

First Name
App Test

Last Name
User

Email
tedis@heartsapiens.art

Phone
+90(555)5555555

Date of Birth
01.01.2000

Gender
☐ Female ☒ Male ☐ Other

Nationality
TR - Türkiye

Target INR Range
☒ 2.00-3.00 ☐ 2.50-3.50

Description

Observation Values

Tedis © 2026 Treatment Monitoring System [Privacy Policy](#)

STEP 5

Notifications Tab

This tab shows in-app notifications sent to you by your doctor or institution (e.g. appointment reminders, general announcements). Notifications are listed with their send time and sender; any attachments can also be viewed here. This is a persistent notification inbox separate from mobile push notifications — if you want to re-read an older notification later, this is where to find it.

**Tediz**
Treatment Monitoring System

Version: 2026.08.01

[Dashboard](#)

[Pending Alerts](#)

[Patient Detail](#)

Gender

☐ Female ☒ Male ☐ Other

Nationality

TR - Türkiye

Target INR Range

☒ 2.00-3.00 ☐ 2.50-3.50

Description

Observation Values

Notifications

Start Time

GG.AA.YYYY

End Time

GG.AA.YYYY

Title	Content	Delivery Status	Send Time	Attachments	Actions
No Records Found					

Tediz © 2026 Treatment Monitoring System

[Privacy Policy](#)

STEP 6

My Account

The "My Account" page shows your login username, full name, email and roles; password fields are always shown empty for security. To change your password, fill in the "Password" and "Confirm Password" fields and click "Update". Further down the page there is also an option to change your email address; this requires entering your current password for security.

Tedis
Treatment Monitoring System
Version: 2026.08.01

[Dashboard](#)
[Pending Alerts](#)
[Patient Detail](#)

My Account [App Test](#) [EN](#) [\[→\]](#)

Account Details

Username
tedis@heartsapiens.art

Password
Password

Confirm Password
Confirm Password

First Name
App Test

Last Name
User

E-Mail
tedis@heartsapiens.art

Roles
☒ Patient
[Edit](#) [Clear](#)

Change Email Address

Current Password

Tedis © 2026 Treatment Monitoring System [Privacy Policy](#)